

**SIMPONI® (subcutaneous golimumab)**  
**Density of SIMPONI Single-Dose Prefilled Syringe and Single-Dose Prefilled SmartJect® Autoinjector**

**SUMMARY**

- The density for SIMPONI 100 mg/mL or 50 mg/0.5 mL single-dose prefilled syringe and single-dose prefilled SmartJect® autoinjector is 1.040 g/mL.<sup>1</sup>

**REFERENCES**

1. Data on File. CNTO148 (golimumab) Module 3.2.P.2.2.3 Drug Product Physicochemical and Biological Properties. Janssen Research and Development; 2008.